

ODDS - Standard License - Standard License

- Central user administration
- DSS Pro recording
- Using meta information like Author IDs for automatic processing
- Email transfer
- Olympus Dictation Portal
- Secure transfer protocols
- Defining recipients for users
- Optional connection to a self hosted ftp-server



Olympus Dictation Delivery Service

The ODDS Standard License grants full access to the Olympus Dictation Delivery Service (ODDS) for one user (12 months). This license allows a user of the Olympus Dictation app to connect it to ODDS and to make use of additional professional dictation features.

The app settings can be centrally administrated via the Olympus Dictation Portal (ODP).

Specifications

Olympus Dictation Portal (ODP)

Central administration	<ul style="list-style-type: none"> Users Licenses Workflow setting Transfer methods Security settings Recipients Keywords ("worktype lists") Recording format Allocating meta information
Languages	<ul style="list-style-type: none"> English German Czech Swedish Spanish French

Olympus Dictation App add-ons

	Olympus Dictation Delivery Service (ODDS) offers advanced workflow and routing options beyond the basic functions of Olympus Dictation App
Automatic processing	Share and send large dictations (files)
Recording Formats	DSS, DSS Pro
DSS Pro file encryption	256 bit AES
Central administration	via ODP
Secure transfer protocols	HTTPS / FTPS / SFTP

Licensing

Duration	12months *
User per license	1
Assignable devices	1 per user
Cancellation	Automatic after expiration
Renewal	Manually

* License period starts after assignment to end-customer account.